



# SNUGGLER

Code no: 06/1270

**“Every farrowing house should have one”**

## **ROTECH SNUGGLER - revive hypothermic piglets**

*The first hours after birth are the most dangerous for piglets, and many suffer stress or even death because of hypothermia. Once they become lethargic their chances of survival severely diminish.*

- **Revive hypothermic piglets**
- **Reduce pre-weaning losses**
- **Easy-to-clean design**

The **Rotech SNUGGLER** comprises a large, high quality insulated box into which a heater module, with digital controller, and a slatted floor have been engineered. Temperature controlled warm air emerges from underneath the slatted floor, ensuring a constant gentle warmth. The box will hold 2 or 3 newborns at a time. In trials, piglets could be removed within an hour.

The **Rotech SNUGGLER** is mains powered, but can be run from a vehicle's 12v supply using a low-cost inverter.

The **Rotech SNUGGLER** has been specially designed with hygiene and bio-security in mind. The heater module simply lifts out to enable thorough pressure washing of the box and floor.

### **Specification:**

*Power Consumption: 230v AC, 120 watts max*

*Weight: 9.4 kg approx*

*Dimensions: 74 x 40 x 46 cm*

*Usable Capacity: 35 Litres approx*

### **Rotech Livestock Equipment Limited**

Woodhorn Lane, Oving, Chichester

PO20 2BX, England

Tel: (+44) (0)1243 787115

Fax: (+44) (0)1243 531847

E-mail: [admin@rotechuk.com](mailto:admin@rotechuk.com)

Website: [www.rotechuk.com](http://www.rotechuk.com)

