



LCD DIGITAL MICROSCOPE

LCD DIGITAL MICROSCOPE with CAMERA

This LCD Digital Microscope combines all of the powerful features of a professional grade microscope with an easy-to-use, fully-integrated, high-resolution 3.5" LCD screen, giving users the ability to view, capture and save vibrant, high-resolution images and videos.

The LCD Digital Microscope replaces traditional eyepieces with an LCD Screen for easy and comfortable viewing for the user and to share with others. It also features a built-in digital camera for quick and easy snapshots or short videos, which can be directly transferred to a computer via the USB cable.

With three objective lens magnifications (x4, x10, and x40), x4 digital zoom, six position colour filter wheel, top and bottom LED illumination, 128MB internal storage memory and built-in SMP digital camera, the LCD Digital Microscope is capable of capturing images and short videos up to **1600 times** their original size.

Specification:

Compound (Biological) Microscope

- USB Cable for transferring images to a PC
- 40 to 400 Power - up to 1600 Power with Digital Zoom
- AC Adaptor to power the microscope
- 3.5" (88mm) LCD Screen with x4 Digital Zoom
- Carrying Case included
- Built-in Digital Camera - 2 Mega Pixels
- Weight 1446g
- Top and bottom LED Illumination
- Two Year Limited Warranty
- Mechanical Stage - 3.5" x 3.5" (88mm x 88mm)
- 128MB Internal Storage Memory
- SD Card Slot
- Five Prepared Slides
- Objective Lens - x4, x10, and x40
- Six Position Colour Filter Wheel (Not available with Heated Stage)



**NOW AVAILABLE
WITH
HEATED STAGE!**



Rotech Livestock Equipment Limited

Woodhorn Lane, Oving, Chichester
PO20 2BX, England

Tel: (+44) (0)1243 787115

Fax: (+44) (0)1243 531847

E-mail: admin@rotechuk.com

Website: www.rotechuk.com