



SEMCKECK 2

Code no: 06/0098

SEMCKECK 2 COLORIMETER

Measure sperm density using the latest technology from Rotech



- Specifically designed for pig AI
- Low cost, easy to operate
- Everything you need, all in one kit
- Reduce boar cost per dose
- Standardise sperm density
- Early warning of reduced boar performance

"GETTING MORE FROM YOUR BOAR"

Rotech Livestock Equipment Limited

Woodhorn Lane, Oving, Chichester
PO20 2BX, England

Tel: (+44) (0)1243 787115

Fax: (+44) (0)1243 531847

E-mail: admin@rotechuk.com

Website: www.rotechuk.com





SEMCKECK 2


Colorimetry — used to measure sperm density

Reduce boar cost per dose

The **Rotech Semcheck 2 Colorimeter** quickly and easily gives an assessment of boar sperm density - enabling the optimum number of semen doses to be prepared from a single collection.

Low cost, easy to operate



The **Semcheck 2** comes complete with easy-to-follow, laminated instructions, including a calibration graph and ready reckoner.



Everything you need, all in one kit:

Rotech can supply the **Semcheck 2** in a kit which includes all the accessories needed to measure semen density:

- Rotech Semcheck 2 Colorimeter
- 0.3ml Fixed Volume Pipetter
- 2.7ml Fixed Volume Pipetter
- Disposable Cuvettes
- 1.0ml Pipetter Tips
- 5.0ml Pipetter Tips



Standardise sperm density

Checking each boar collection with the **Semcheck 2 Colorimeter** enables sperm density to be standardised in all semen doses.

Early warning of reduced boar performance

Colorimetry on a regular basis will highlight the earliest signs of any reduced boar performance.

Specifically designed for pig AI

The **Semcheck 2 Colorimeter** has been specifically designed for **pig AI**.

Rotech have many years of experience with AI labs using colorimeters and can provide all the necessary training, instruction and advice.

It should be remembered that determining sperm density using a colorimeter is not 100% accurate and should only be used as a guide.

Rotech Livestock Equipment Limited

Woodhorn Lane, Oving, Chichester
PO20 2BX, England

Tel: (+44) (0)1243 787115

Fax: (+44) (0)1243 531847

E-mail: admin@rotechuk.com

Website: www.rotechuk.com

