



MICROSCOPE SLIDE-WARMER

Code no: 06/0045

WARMS MICROSCOPE SLIDES FOR SEMEN MOTILITY EXAMINATION

An economical alternative to a heated stage, the Rotech Slide Warmer will warm a handful of slides to 35°C for motility examination of semen.

Semen motility examination must be carried out using a warm slide, because low temperature will have a negative effect on motility, leading to the possibility of an inaccurate motility score being given.

A heated stage will keep a slide warm during the examination, but it is not practical for warming it prior to preparing the sample. The Slide Warmer simply needs to be plugged into the mains supply and slides placed on its surface. The thermostatically controlled device will then maintain a temperature of 35°C, so that when you come to place a drop of semen for examination, there will be no temperature shock. Because motility checking is so quick, there may even be no need to have a heated stage at all.

Specification:

Dimensions: 155 x 155 x 25 mm

Power requirements: 220v or 110v AC (please specify when ordering)

Surface temperature: Pre-adjusted to 35°C.



Rotech Livestock Equipment Limited

Woodhorn Lane, Oving, Chichester
PO20 2BX, England

Tel: (+44) (0)1243 787115

Fax: (+44) (0)1243 531847

E-mail: admin@rotechuk.com

Website: www.rotechuk.com

