



THERMOSAFE 25

Code no: 06/1230

MORE THAN JUST A 17°C SEMEN STORAGE BOX:

- SEMEN STORAGE 17°C
- REFRIGERATOR 4°C
- WARM-BOX 35°C

Re-programmable, digitally controlled temperature

Not only will the **Thermosafe 25** maintain boar semen accurately at its optimum storage temperature of **17°C**, **but** it can easily be re-programmed by the user to run at any temperature, up to approximately **35°C** above or **20°C** below ambient.*

Therefore, for example, you could re-set the digital controller to run the **Thermosafe 25** at **4°C** for storing vaccines or other medicines, or to **35°C** for warming collection equipment, microscope slides, or doses of semen to blood temperature.

The **Thermosafe 25** will run on 12 volts, from a vehicle power socket, and it incorporates a transformer for running from the mains.

34 Litre nominal capacity, will hold approx 80-100 doses, depending upon type of pack.

With "PELTIER PROTECTION" to maximise the life of the peltier heating and cooling element.

Accessories:

Battery saver device - ensures that your vehicle battery is protected even if the Thermosafe is left turned on whilst the engine is off.

Storage Baskets - to help separate semen doses and promote improved airflow.

* Tested using a mains-powered transformer. Performance when run from a 12v battery may be lower, depending on battery condition.



Rotech Livestock Equipment Limited

Woodhorn Lane, Oving, Chichester
PO20 2BX, England

Tel: (+44) (0)1243 787115

Fax: (+44) (0)1243 531847

E-mail: admin@rotechuk.com

Website: www.rotechuk.com

