



WARMING BOX

Code no: 06/1260

Warm, or dry, lab equipment with the:

Rotech WARMING BOX

An essential part of pig AI is the protection of raw semen from temperature shock. Everything which may come into contact with raw semen should therefore be heated to blood temperature, 37°C.



- **Simple to use**
- **Low cost**
- **Digitally controlled, even temperature throughout**
- **Clear LED internal temperature display**
- **Mains powered**

Rotech pioneered the move towards disposable, single-use items for pig AI. Previously everything had to be sterilised and then dried in large warming cabinets. Nowadays, using disposables means no more sterilising, but small quantities of equipment (collection flasks, thermometers, microscope slides etc) must be warmed, and then kept warm, for collection purposes.

Specification:

Power Consumption: 230v AC, 70 watts max

Weight: 6.5 kg approx

Dimensions: 66 x 37 x 37 cm approx

Usable Capacity: 24 Litres approx



Rotech Livestock Equipment Limited

Woodhorn Lane, Oving, Chichester

PO20 2BX, England

Tel: (+44) (0)1243 787115

Fax: (+44) (0)1243 531847

E-mail: admin@rotechuk.com

Website: www.rotechuk.com

