



# WARMING MATS

*Keeps diluent at 35°C on your worktop*

## HEATED WORK SURFACES KEEP DILUENT WARM

A **Rotech WARMING MAT** will transform your work surface into a temperature-controlled environment. It will keep your diluent at 35°C in the Sembag™ .

- **Constant, easy-to-read LED temperature display**
- **Digitally controlled temperature setting, simple to use**
- **Warm several Sembags™ at once**
- **High quality, wipe-clean surface**
- **Low cost, convenient + efficient**
- **Range of sizes to suit requirements**



Can also be used to warm microscope slides, cover slips, thermometers etc

*Requires 220v AC mains supply.*

### Product Range:

06/1185	Warming Mat	Size 300mm x 900mm
06/1186	Warming Mat	Size 300mm x 1600mm
06/1187	Warming Mat	Size 600mm x 1600mm

*Also available made to order to suit specific size requirements*

### Rotech Livestock Equipment Limited

Woodhorn Lane, Oving, Chichester  
PO20 2BX, England

Tel: (+44) (0)1243 787115

Fax: (+44) (0)1243 531847

E-mail: [admin@rotechuk.com](mailto:admin@rotechuk.com)

Website: [www.rotechuk.com](http://www.rotechuk.com)

