



# HEATED WATERBATH

Code no: 06/1166

The **Rotech Heated Waterbath** has been especially designed with the needs of the small farm **Artificial Insemination** laboratory in mind.

The waterbath uses a digital controller to ensure that the water temperature remains constantly at the chosen temperature, normally 35°C, more or less the blood temperature of the animal. The controller can easily be re-set to operate at a different temperature, up to 60°C, if required for other purposes. However since it is of plastic construction, it will not support boiling water.



#### Usable internal dimensions -

*above perforated floor -*  
36 cms x 26 cms x 18 cms  
(16.5 Litres approx)

#### Specification -

Power requirements: 230V  
Max operating temp: 60°C  
Temperature re-set using  
push-button digital controller

Using an electric kettle element and a digital controller, the **Rotech Waterbath** is a simple, robust and accurate piece of equipment which will give a performance far beyond what might be concluded from its price tag.

Each bath is carefully calibrated and tested to 35°C during manufacture. Re-setting to run at different temperatures is quick and easy, and full instructions are provided.

#### Rotech Livestock Equipment Limited

Woodhorn Lane, Oving, Chichester  
PO20 2BX, England

Tel: (+44) (0)1243 787115

Fax: (+44) (0)1243 531847

E-mail: [admin@rotechuk.com](mailto:admin@rotechuk.com)

Website: [www.rotechuk.com](http://www.rotechuk.com)

